

Order No.: DD+DIS072.02E



1 piece UXWJU MA1

ADC Compact Plus

Type 5146 / 100/200

Rating			When?		
safety relevant	reliability	new function	immediately	next service	as required
		X			X

Please file this document in section 10 of your technical documentation

Device Software ACP_1010 available on MedNet

Device Software ACP_1010 can be downloaded from MedNet as follows:

- (1) Open your browser and go to the MedNet CSO Library
- (2) Select
 Computed Radiology
 ADC Compact Plus
 CSO Freeware
- (3) Select
 Device Software ACP1010
- (4) Download the attachment file Device Software ACP1010.zip
- (5) Unzip Device Software ACP1010.zip
- (6) Copy the following directories to one floppy each and label them as described:

Disk name:	Label name:
disk1	→ ACP1010_1_3
disk2	→ ACP1010_2_3
disk3	→ ACP1010_3_3
harddisk formatter	→ 4
Portex 1	→ 5
Portex 2	→ 6



How to label a floppy disk:

1. Insert floppy and select floppy drive in Explorer
2. Click right mouse button and select Properties
3. Select tab General and enter required label name